

LabArchives is invested in supporting digital learning and Open Educational Resources (OER). Hundreds of labs and protocols are available through Lab Builder, our open source library.

To find out more and to gain free access to our complete library of content, please visit our [website](#) to sign up for an account.


Lab Builder, by LabArchives, also offers OpenStax textbooks at no cost inside of your LabArchives account. Build a complete course using the LabArchives notebook platform, course management tools, and our open digital library. Please contact [Support](#) with any questions on how begin implementation in to your course.



**Find precisely what you need in our digital library for your course.**



**Find, customize, and build content for your course using the LabArchives notebook platform and its course management tools.**



To stay up to date on all things LabArchives and read about our numerous success stories, please visit our [Blog](#).



## Agarose Gel Electrophoresis Protocol

Dan Kraut, Ph.D.  
Department of Chemistry and  
Biochemistry  
Villanova University

Julie Wolf, M.S.  
Department of Biological Sciences  
University of Maryland Baltimore County

Natalie Stringer, Ph.D.  
Education Division  
LabArchives, LLC

### Attribution

Molecular Biology Course Pack by LabArchives is licensed under a [Creative Commons Attribution-NonCommercial-ShareAlike 4.0 International License](https://creativecommons.org/licenses/by-nc-sa/4.0/).

You must give appropriate credit, provide a link to the license, and indicate if changes were made. You may do so in any reasonable manner, but not in any way that suggests the licensor endorses you or your use. All content included in Lab Builder includes rights to use within LabArchives. For detailed information on various Creative Commons type licenses for the content which appear in Lab Builder go to [www.creativecommons.org](http://www.creativecommons.org).



### Topic

- PCR
- Electrophoresis
- DNA
- Mutations
- Genetics
- Agarose
- Molecular weight standards

### Timeline

- **Prepare samples, run and image gel:** 45 minutes – 1.5 hours

### Safety

- Wear gloves and safety goggles!

## Protocol – Molecular Biology: Agarose Gel Electrophoresis

### Equipment

- Gel casting tray
- Combs (to make the sample wells)
- Gel box to hold buffer and gel
- Voltage source
- Imager (many options)
- Microwave

### Materials

- TAE or TBE
- Agarose
- [Molecular weight standards](#)
- DNA staining reagent
  - Ethidium bromide
  - SYBR Safe
  - GelRed or GelGreen





Image credit: Dan Kraut, Villanova University

### Prepare the Agarose Gel

1. Prepare casting tray
2. Prepare 0.7-1.5% agarose solution in 1X TBE or TAE buffer
  - Use the same buffer to run your gel
  - The concentration depends on the size of the DNA you need to separate
  - The final volume of agarose depends on your specific gel apparatus
3. Mix agarose with buffer in a microwavable flask



## Protocol – Molecular Biology: Agarose Gel Electrophoresis

4. Microwave until the solution bubbles (without overflowing) and the agarose is completely dissolved (about a minute is a good place to start).
5. Let solution cool
  - Be sure to use paper towels or a hot pad to avoid burning yourself
  - Allow to cool until you can touch the flask with a gloved hand without it being too hot to handle; you can use a cold water bath or run the flask under tap-water to help.
6. Add DNA staining reagent to appropriate final concentration and swirl
  - *Alternative* - stain after running the gel
7. Pour the solution into the casting tray
8. Allow to solidify

### Running the Gel

1. Add loading buffer to each of your samples in an eppendorf tube
2. Amount of DNA loaded on gel depends on many [factors](#)
  - Fragment size
  - Detection method
  - Well volume
  - Distribution of bands (fragment sizes)
3. Add the appropriate amount of sample to each well of your gel
4. Add the vendor's recommended volume of molecular weight standards (ie Ladder)
5. [Run gel](#)
  - Optimal voltage and time depend on fragment size and buffer choice